

Rejuvenate™ 2.0 Hybrid Copper Mattress Presentation



www.nscmedical.com



Rejuvenate™ 2.0 Hybrid Copper



Elevate Your Sleep with the Rejuvenate™ 2.0 Hybrid Copper Series

Introduction

Experience the transformative sleep on our Rejuvenate™ 2.0 Hybrid Copper series, featuring innovative self-sterilizing technology. This consumer sleep surface incorporates a patented copper contact layer, recognized as a Class 1 medical device for biological change, setting new standards for health-conscious sleep environments.



Rejuvenate™ 2.0 Hybrid Copper




- 99% bacteria kill rate within 24 hours through your mattress protector, bed sheets, and sleeping apparel.
- Apparel assists with inflammation reduction and virus reduction.

Layers Composition

1. Copper Contact Surface:

- Manufactured in a GOTS certified facility in Canada.
- Recognized as a Class 1 medical device for biological change.



Rejuvenate™ 2.0 Hybrid Copper



2. Natural Latex Layer:

- Hypoallergenic and all-natural.
- ABC Laboratories for antimicrobial.
- GOLS certified for 100% natural content.
- ECO institute certified for 100% natural latex.
- STANDARD 100 by OEKO-TEX® certified for natural content.


3. Medi-Gel Layer:

- Assists with spinal decompression, pressure point relief, and motion transfer.
- Clinical Medi-Gel Technology.
- Outperforms other air and gel technologies based on pressure mapping.



4. Medi-Foam Base:

- Assists with weight-bearing support in all sleeping positions.
- NSC Medical's exclusive formulation neutralizes off-gassing smell and reduces heat retention.
- Soy-based alternatives reduce carbon footprint by up to 20% and are Certi-Pur certified.



Rejuvenate™ 2.0 Hybrid Copper



Features

- Low Volatile Organic Compound (VOC) soy emissions for indoor air quality (less than 0.5 parts per million).
- Removable Waterfall copper cover with piping around the top perimeter.
- Handmade mattress with a manufacturing variance of 2.5 cm (1 in.) on the crown height.
- Recommended for use with posture box, box spring, slat system (with space no more than 3 inches), or adjustable bed.
- Allow 14 to 21 days for NSC Medical mattresses to adapt to your body.

Developed and manufactured in Canada.



Certificates

- CertiPUR-US® certified soy foam: Free from ozone depleters, flame retardants, heavy metals, formaldehyde, and phthalates.
- ISO-9001:2015 Accredited Manufacturing Processes.
- ISO-14001:2015 Accredited Manufacturing Processes.
- Manufacturer is Class I medical device certified.
- Manufacturing traceability on all components ISO-9001:2017 accredited.

Rejuvenate™ 2.0 Hybrid Copper



In Conclusion

The Rejuvenate™ 2.0 Hybrid Copper mattress is more than a sleep surface; it's a commitment to your well-being. Elevate your sleep with this innovative hybrid, where technology meets comfort. Invest in your sleep; invest in your health.



Check out our mattress models on our website by clicking the link below. You'll find purchase links there too!

www.nscmedical.com



www.nscmedical.com